VELTEK ASSOCIATES, INC.

TECHNICAL DATA FILES

ALCOH-WIPE®

70% USP Isopropyl Alcohol and 30% Water for
Injection

Individually Packaged and Pre-Saturated 100%
Polyester Wipe

29SEP2022

Copyright © 2022 Veltek Associates, Inc.

All Rights Reserved. Confidential and Proprietary Property of Veltek Associates, Inc. No part of this document may be reproduced or disclosed, in whole or in part, without the specific written permission of Veltek Associates, Inc.

TABLE OF CONTENTS

PRODUCT OVERVIEW.....3
 Product Uses3
 Quality and Manufacturing4
 Features and Benefits4
CHEMICAL SPECIFICATIONS5
MATERIAL SPECIFICATIONS5
ORDERING INFORMATION5
LABELING6
ADDITIONAL DOCUMENTATION8

ALCHO-WIPE VEL6-12X12-S-2302

PRODUCT OVERVIEW

The **ALCOH-WIPE** is ready-to-use, individually packaged, and saturated with VAI®'s 70% USP Isopropyl Alcohol (IPA) and 30% USP Water for Injection. These wipes can be used in both aseptic and non-aseptic wipe downs of filling and packaging machinery, stainless steel surfaces, lexan, polycarbonate, glass, and any critical surface that requires cleaning. The **ALCOH-WIPE** has been specifically designed for operations that require the use of an individually packaged and sterile IPA wipe. Pre-saturated wipers, including **ALCOH-WIPE**, can prevent the overuse of chemicals and reduce costs. The **ALCOH-WIPE** is available sterile in 6" × 6" (15.2 cm × 15.2 cm), 12" × 12" (30.5 cm × 30.5 cm), and 18" × 18" (45.7 cm × 45.7 cm) sizes.

ALCOH-WIPE is filled in ISO 5 certified cleanrooms to ISO 7 operations (Grade A/B, Former Class 10,000), filtered at 0.2 microns, and subsequently terminally sterilized to 10⁻⁶ sterility assurance level. Individual lots are lot sterility tested according to current USP compendium, are completely traceable, and completely validated for sterility and shelf life. **ALCOH-WIPE** is delivered each time with lot specific Certificate of Analysis, Certificate of Irradiation, Certificate of Sterility, and Certificate of Conformance tested to current USP compendium.

The material used in our **ALCOH-WIPE** is a 100% continuous knit polyester filament that has been specifically designed to be exceptionally clean and absorbent. Particulates, residues, and foreign matter is effectively entrapped in the webbing of the material. The knit construction and pattern has been developed to ensure the greatest product integrity and cleaning performance while offering high resistance to abrasion and picking.

FocusEdge cutting technology is the standard that VAI uses to cut their wipes. This cutting method does not heat the yarn bundles as high as lasers, therefore, eliminating the small plastic balls that form on the end of every yarn strand. Due to the wipe's sealed edges, shedding is minimal during wipe downs or general use, therefore, reducing particulates in controlled environments. After cutting, the wipes are then laundered in ISO Class 4 conditions to ensure cleanliness.

VAI's **ALCOH-WIPE** is packaged in an opaque, ported, and easy-tear style bag. Each wiper is individually bagged and packed in two liner bags for easy introduction into the core or other critical environments.

Product Uses


- Aseptic or non-aseptic environments
- Wipe downs on filling and packing machinery
- Spill control
- On lexan, glass, polycarbonate, stainless steel, or any other critical surface that requires cleaning
- In operations that require an individually bagged IPA wiper
- To prevent overuse of IPA when compared to using trigger sprayers or gallon pour bottles

Quality and Manufacturing

- Wipe material laundered in ISO 4 (Former Class 10)
- 70% USP Isopropyl Alcohol and 30% USP Water for Injection (<0.25 EU/mL) solution is filtered at 0.2 microns
- Filled in ISO 5 certified cleanroom to ISO 7 operations (Grade A/B, Former Class 10,000)
- Cut using FocusEdge cutting technology
- Gamma irradiated to 10^{-6} SAL
- Lot sterility tested according to current USP Compendium
- Completely lot traceable
- Completely validated for sterility and shelf life

Features and Benefits

- ISO 5 certified cleanroom to ISO 7 operations (Grade A/B, Former Class 10,000) Cleanliness
- Individually bagged with two liner bags
- Packaged in an opaque, easy-tear style bag
- Single wipe is folded for easy removal
- Individually labeled with lot number and expiration
- For use on multiple surfaces
- Low endotoxins
- Reduces VOC emissions
- 100% synthetic filament polyester constructed in a continuous knit that minimizes fiber and particulate release
- Material is exceptionally clean and absorbent
- Resistance to abrasion and picking
- FocusEdge allows for minimal shedding via sealed edges
- Ready-to-use
- Delivered with lot specific Certificate of Irradiation, Certificate of Sterility, Certificate of Analysis, and Certificate of Conformance tested to current USP compendium

100% Polyester Wiper Substrate

CHEMICAL SPECIFICATIONS

ALCOH-WIPE – 70% USP Isopropyl Alcohol and 30% Water for Injection	
Appearance	Clear Colorless
Assay	68.0 – 72.0%
Acidity	<= 1.0 mL
Nonvolatile Residue	<5.0 mg
Specific Gravity @ 20 Degrees C	0.872 – 0.883
Expiration Period	2 years

MATERIAL SPECIFICATIONS

ALCOH-WIPE – Laundered Polyester Wipe	
Material Description	Laundered Polyester Wipe
Dimensions	31.30 cm × 31.50 cm
Sorbancy	>= 270 mL/m ²
Particles	<4.5 × 10 ⁶ particles/m ²
Material	100% Polyester
Construction	Continuous Filament Knitted Construction
Color	White
Classification Compatibility	ISO 5 certified cleanroom to ISO 7 operations (Grade A/B, Former Class 10,000)
Thickness	0.56 – 0.64 mm
Fibers @ 100 um/m ²	63, 1.0 IEST-RP-CC004.3, SEC 6.1.3/6.2.2

ORDERING INFORMATION

ALCOH-WIPE – 70% USP Isopropyl Alcohol and 30% Water for Injection		
Part Number	Description	Qty/Case
VEL6-6X6-S-2307	ALCOH-WIPE, 6" × 6" (15.2 cm × 15.2 cm), Individually Bagged, Sterile	100
VEL6-12X12-S-2302	ALCOH-WIPE, 12" × 12" (30.5 cm × 30.5 cm), Individually Bagged, Sterile	100
VEL6-18X18-S-2304	ALCOH-WIPE, 18" × 18" (45.7 cm × 45.7 cm), Individually Bagged, Sterile	100

Veltek Associates, Inc.

15 Lee Boulevard, Malvern, PA 19355-1234
 Tel: 610-644-8335 Fax: 610-644-8335 www.sterile.com
 29SEP2022

LABELING

ALCOH-WIPE®
70% USP Isopropyl Alcohol and 30% Water for Injection

DANGER GEFAHR PERIGO PELIGRO
PERICOLO 危險 FARE FARA 위험

EN
70% Isopropyl Alcohol USP formulated with USP / EP WFI grade water

INGREDIENTS:
 Isopropanol (CAS 67-63-0) 70%
 Other ingredients: water

DANGER. Highly flammable liquid and vapour. Causes serious eye irritation. May cause drowsiness or dizziness.

PRECAUTIONARY STATEMENTS: Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. Use only outdoors or in a well-ventilated area. Wear eye protection, protective gloves, protective clothing. In case of fire: Use alcohol resistant foam, carbon dioxide (CO₂), Dry powder to extinguish. Store in a well-ventilated place. Keep cool. Dispose of contents and container to an authorised waste collection point.

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6 in x 6 in (15,25 cm x 15,25 cm) 18 in x 18 in (45,75 cm x 45,75 cm)
 12 in x 12 in (30,5 cm x 30,5 cm)

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Red dot indicates Sterile
 See the container for the lot #/ expiration date

SDS: VEL-038-AG (Americas), VEL-038-AG-EU (Europe)
Emergency Number: 1-866-928-0789 (toll free USA), 1-800-579-7421 (toll free Canada), 1-215-207-0061 (Americas), +44 1235 239670 (Europe), +44 1235 239671 (Middle East/Africa) +353 (0) 809 280 (Ireland, available to the public: 08:00-22:00), 000 800 100 7479 (Hind).

See Safety Data Sheet for first aid and additional product information.

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KO (한국어)
ES | NL | SV

Case Label

Any specific product label is available upon request.

ALCOH-WIPE®

70% Isopropyl Alcohol USP formulated with USP/EP WFI grade water

DANGER

INGREDIENTS

Isopropanol (CAS 67-63-0)..... 70%

Other Ingredients: water

DANGER

Highly flammable liquid and vapour. Causes serious eye irritation. May cause drowsiness or dizziness.

PRECAUTIONARY STATEMENTS

Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. Use only outdoors or in a well-ventilated area. Wear eye protection, protective gloves, protective clothing. In case of fire: Use alcohol resistant foam, carbon dioxide (CO₂), Dry powder to extinguish. Store in a well-ventilated place. Keep cool. Dispose of contents and container to an authorized waste collection point.

SDS: VEL-038-AW

Emergency Number

1-866-928-0789 (toll free USA), 1-800-579-7421 (toll free Canada), 1-215-207-0061 (Americas), +52-55-5004-8763 (Mexico).

See Safety Data Sheet for first aid and additional product information.

For Industrial Use Only

Processed, Filled, Packaged, Sterilized, and Tested in the USA.

Dry Wipes Made in South Korea.

Patents: www.sterile.com/patents

Manufactured By

Veltek Associates, Inc.

15 Lee Blvd.

Malvern, PA 19355-1234 USA

Tel: 610-644-8335

www.sterile.com

Veltek Associates, Inc.

15 Lee Boulevard, Malvern, PA 19355-1234
Tel: 610-644-8335 Fax: 610-644-8335 www.sterile.com
29SEP2022

ADDITIONAL DOCUMENTATION

Upon request, the following additional documentation is available:

- Product Testing Reports
- Product Labels
- Product Validation
- Safety Data Sheet
- Sample Lot Specific Documentation – Product Release Certificate that includes: Certificate of Sterility, Certificate of Irradiation, Certificate of Analysis, and Certificate of Conformance

VAI®'s Wiper Family

Veltek Associates, Inc. *Disposable Products Manufacturing Division's* (DPMD) capabilities include sterile disposable garments, sterile disposable face masks, dry wipers, sterile wipers, custom cleanroom documentation, sterile synthetic writing substrate, and a HEPA filtered cleanroom printer. DPMD has responded and met demands of professionals in the industry since 1985. Our manufacturing capabilities assure that we will be able to meet all requests and fulfill all of our customer's needs. All DPMD products are processed and packaged in cleanroom facilities from beginning to end and subsequently gamma irradiated or aseptically filled in order assure sterility. For more information call 610-644-8335, email sales@sterile.com or visit our website at www.sterile.com. Patents: www.sterile.com/patents

Veltek Associates, Inc.

15 Lee Boulevard, Malvern, PA 19355-1234
Tel: 610-644-8335 Fax: 610-644-8335 www.sterile.com
29SEP2022