



PendoTECH® Single Use Pressure Sensors™ for Use in Flexible BioProcess Containers

- Easily installs onto a flexible bioprocess container
- Measures pressure in the container
- Accurate performance in either liquid or gas applications
- No calibration required
- Low cost for single use applications
- Qualified for use up to 15psi (1.03bar)
- Precision better than ± 0.003 psi (0.2mbar) with PressureMAT® High Resolution Monitor



Applications Include:

- Safety: Preventing disposable bioreactor and storage bag over-pressurization
- Performance: Disposable bioreactor monitoring
- Leak Testing: For disposable bag pre-use integrity testing[^]
- Volume Measurement: Measuring liquid level by hydrostatic pressure (Pascal's Principle)[^]

Port Plate Details:

- Port plate with sensor receptacle seals to the film during the manufacturing process of the container
- Made of polyethylene specifically for port plate use
- Compatible with gamma irradiation
- NaOH resistant
- Meets USP Class VI

Sensor Details:

- 100% tested for accuracy and leaks during manufacturing
- Available in polycarbonate or caustic resistant polysulfone material of construction
- All materials meet USP Class VI requirements
- Certificate of quality included with lot certification; individual NIST Calibration Certificates are optional
- Manufactured in an ISO Class 7 clean room
- Compatible with gamma irradiation
- May be non-invasively tested in-place via test port


Product Details:

- After manufacture of the container, the sensor is inserted into the port plate receptacle.
- Locking collar holds sensor securely in place
- Features a double o-ring seal & a stop to ensure the sensor is fully inserted




Final Assembly of
Sensor and Port Plate

Two o-rings to create
seal within the receptacle




Sensor

Stop to indicate
sensor is fully inserted




Sensor Test Port




[^] Link to article: PendoTECH Single Use Pressure Sensors: Calibration, Accuracy, and Implementation
PendoTECH® and PressureMAT® are registered trademarks of PendoTECH, all rights reserved.

Selection of Monitors to Read Sensors



PressureMAT® Monitor, Alarm & Transmitter



PressureMAT® Sensor Transmitter



PendoTECH PressureChecker™ Non-invasive Sensor Verification Tool

Installation Example



Sensor connects to monitor for reading

Tote with container installed with sensor cable lead exiting through hole in bottom of tote

Zoom view of sensor affixed to port plate at bottom of container



Sensor Specifications:

Detail	Specifications																
Manufacturing Testing	100% of sensors are tested for critical quality attributes. - Each sensor is leak tested on the liquid side to confirm integral assembly - Sensors with a test port are leak tested on the test port side to confirm proper atmospheric reference - Each sensor is tested electrically to confirm proper electrical performance - Each sensor is tested to be accurate at 60 psi (4.14 bar) within +/-5% of reading (+/-3.0 psi/0.21bar)																
Accuracy	<table border="1"> <thead> <tr> <th>Positive Range</th> <th>Specification</th> <th>Vacuum Range</th> <th>Specification</th> </tr> </thead> <tbody> <tr> <td>0 to 6psi</td> <td>±2% of reading</td> <td>0 to -7psi</td> <td>±3% of reading</td> </tr> <tr> <td>6 to 30psi</td> <td>±3% of reading</td> <td>-7 to -10psi</td> <td>±5% of reading</td> </tr> <tr> <td>30 to 60psi</td> <td>±5% of reading</td> <td></td> <td></td> </tr> </tbody> </table>	Positive Range	Specification	Vacuum Range	Specification	0 to 6psi	±2% of reading	0 to -7psi	±3% of reading	6 to 30psi	±3% of reading	-7 to -10psi	±5% of reading	30 to 60psi	±5% of reading		
Positive Range	Specification	Vacuum Range	Specification														
0 to 6psi	±2% of reading	0 to -7psi	±3% of reading														
6 to 30psi	±3% of reading	-7 to -10psi	±5% of reading														
30 to 60psi	±5% of reading																
Pressure range	For full assembly: up to 15psi For insert only: 75psi																
Biocompatibility	All materials in contact with product fluid path meet USP Class VI requirements, both pre and post irradiation																
Manufacturing Environment	ISO 9001 certified facility; Class 7 clean room																
Gamma Irradiation	Up to 50 kiloGrays																
X-ray Irradiation	Up to 50 kiloGrays																
Operating temperature	2°C to 40°C (other ranges with process qualification)																
Storage temperature	-25°C to 65°C																
Input/Output impedance	270 Ohms to 400 Ohms																
Excitation voltage	2.5 to 10 volts DC (for best long term stability, use a lower excitation voltage)																
Sensor Output	0.2584 mV/Volt/psi																
Connector	Rating: IP67 when connected to reusable cable																
Shelf life	5 years																
Packaging	White Tyvek® and clear pouch with easy-open chevron seal; box of 25 sensors in polyethylene bags																

Ordering Information

PTPL-PRESS	Single Use Pressure Sensor Insert for Port Plate, Polycarbonate
PTPL-PREPS	Single Use Pressure Sensor Insert for Port Plate, Polysulfone
PORT-D1-HDPE	Port Plate for Sensor Inserts, 2.25"D, Material: HDPE DOW DMDA-8007 HEALTH+. Pack of 25
PORT-RING	Pack of 25 locking rings for PORT-D1-HDPE bagged in polyethylene bag
PDKT-650-950	PressureChecker pressure sensor and monitor verification tool, psi
Sensors with NIST calibration certificates, see primary pressure sensor datasheet for further details	
PTPL-PREPS-B	Single Use Pressure Sensor Insert for Port Plate, non-sterile, Polysulfone, with 0-6 psi NIST certificate
PTPL-PRESS-B	Single Use Pressure Sensor Insert for Port Plate, non-sterile, Polycarbonate, with 0-6 psi NIST certificate
PTPL-PREPS-W	Single Use Pressure Sensor Insert for Port Plate, non-sterile, Polysulfone with IP67 Connector

Covered by Pat. 10,215,598; Pat. 10,041,896; Pat. EP 14867894.9; and Pat. Pending
 For warranty information see our website at <http://www.pendotech.com/warranty>