

EZwaste[®] Heavy Duty **Solvent Waste System**

www.foxxlifesciences.com



ISO 13485:2016



EZwaste[®] Heavy Duty Solvent Waste System

www.foxxlifesciences.com


EZwaste Solvent Waste System

0			and the series that the	
	1/8" OD PORTS	1/4" OD PORTS	1/4" or 3/8" HB PORT	PART NUMBER
2.5L	6	-	-	332-1412-OEM
	6	-	1	332-1418-OEM
	4	3	-	332-1413-OEM
	4	4	-	332-1421-OEM
	4	3	1	332-1419-OEM
4L	6	-	-	332-9412-OEM
	6	-	1	332-9418-OEM
	4	3	-	332-9413-OEM
	4	4	-	332-9421-OEM
	4	3	1	332-9419-OEM
2T	6	-	-	332-2412-OEM
	6	-	1	332-2418-OEM
	4	3	-	332-2413-OEM
	4	4	-	332-2421-OEM
	4	3	1	332-2419-OEM
10L	6	-	-	332-3412-OEM
	6	-	1	332-3418-OEM
	4	3	-	332-3413-OEM
	4	4	-	332-3421-OEM
	4	3	1	332-3419-OEM
20L	6	-	-	332-4412-OEM
	6	-	1	332-4418-OEM
	4	3	-	332-4413-OEM
	4	4	-	332-4421-OEM
	4	3	1	332-4419-OEM
20L ROUND	6	-	-	332-D412-OEM
	6	-	1	332-D418-OEM
	4	3	-	332-D413-OEM
	4	4	-	332-D421-OEM
	4	3	1	332-D419-OEM


* 1/16" OD Adapters are also avaliable for the HD line as cap only assemblies.

Solvent Waste System Sizing Guide

Due to all the possibilities with configuring your EZwaste Solvent Waste System, we have provided the following images to help identify the fitting required for your laboratory setup. All measurements are approximate.



For measuring ID, place end of tube flat on the piece of paper and determine if the dark area is contained by the tubing.



Rev. 33-001-03 05/23/2016 Foxx Life Sciences, 6 Delaware Drive, Salem, NH 03079 www.foxxlifesciences.com Order Now: 1 (603) 890-FOXX (3699)

