


# VersaCap®, VersaBarb®, and Spigot Technology foxxlifesciences.com



- VersaCap® technology provides the flexibility of interchangeable top-connectors
- VersaCap® technology is offered in 38-430mm, 53B, 83B, and 120mm industry standard sizes
- Quick-connect adapters available for easy port changes
- · Leak-proof cap with replaceable gasket
- VersaCap® design enables stacking
- Available in a variety of configurations
- VersaCap<sup>®</sup>, VersaBarb<sup>®</sup>, and Spigots are autoclavable
- VersaBarb® offers a convenient solution for bottom-connections and can fit a variety of hose sizes
- Spigot offers two flow rates

Foxx Life Sciences has enhanced the features and benefits of the **EZ**grip<sup>™</sup> carboys and bottles by offering the **Versa**Cap<sup>®</sup>, **Versa**Barb<sup>®</sup>, and state-of-the-art spigot technology. These three components expand the wide array of applications supported by the **EZ**grip<sup>™</sup> carboys and bottles.

Patented VersaCap® is leak-proof, has replaceable platinum-cured silicone gaskets, and provides the ultimate in top-connection flexibility. Open VersaCap® is compatible with a variety of adapter inserts, enabling users to quickly change the top connections of each carboy or bottle to match their application needs. VersaCap® offers large handles with finger grips that allow for easy tightening and opening. They are also stackable and are offered in 38-430mm, 53B, 83B (compatible with most carbovs), and 120mm wide-mouth versions that fit exclusively on Foxx Life Sciences' 20L, 40L, and 75L EZFil® carboys.

Quick-connect hose barb options include quick-connect port combinations of two, three, and four, with 1/4" and 3/8" hose barb sizes, as well as molded-in hose barb options, available with  $\frac{1}{8}$ ",  $\frac{1}{4}$ ",  $\frac{5}{16}$ ", and  $\frac{1}{2}$ " sizes. The Clearview adapter is also available for users who need a crystal-clear window into the carboy for visualizing solution homogenization.

For bottom tubulation options, VersaBarb® provides flexibility in drainage, pump connection, filling bell, and automated equipment. VersaBarb® is available in 1/4", 3/8", 1/2", and 3/4" hose barb sizes. A 3/4" sanitary fitting **Versa**Barb® is also available to enable the EZgrip™ carboy to function as a supply reservoir to a larger system.

Foxx Life Sciences' **EZ**grip<sup>™</sup> carboy spigots have large levers that control two flow rates, one for continuous. locked flow and the other for a faster, manual flow. The spigot spout features 5/8" and 3/4" hose barbs, and a standard internal luer tape to fit a male barb adapter. The spigot body has internal threads with a hard stop to prevent over-tightening which may damage the seal.

All VersaCap®, adapters, VersaBarb®, and spigots are made from polypropylene, which may be autoclaved at 15 PSI and 121°C (250°F) for 20 minutes.



### VersaCap® Technology

foxxlifesciences.com



#### VersaCap® 38-430mm



# 205-0123-RLS Cap, Closed Cap Size No. of Quick Connects Qty/pk 38-430mm n/a n/a 4



205-0124	4-RLS			
Cap, Open	with closed ac	dapter		
Cap Size	Hose Barb Size	No. of Quick Connects	Qty/pk	
38-430mm	n/a	n/a	4	
	Cap, Open Cap Size	Cap Size Hose Barb Size	Cap, Open with closed adapter Cap Size Hose Barb No. of Quick Size Connects	Cap, Open with closed adapter Cap Size Hose Barb No. of Quick Qty/pk Size Connects



205-0129-RLS			
Cap, Open	top with adap	ter, ¼" HB x 2	
Cap Size	Hose Barb Size	No. of Quick Connects	Qty/pk
38-430mm	1/4"	n/a	4



205-0120	D-KLO		
Cap, Open	top with adap	oter, 1/8" HB x 3	
Cap Size	Hose Barb Size	No. of Quick Connects	Qty/pk
38-430mm	1/8"	n/a	4













### VersaCap® Technology

foxxlifesciences.com



#### VersaCap® 53B



205-0313-RLS				
Cap, Close	ed			
Cap Size	Hose Barb Size	No. of Quick Connects	Qty/pk	
53B	n/a	n/a	2	











205-333	1-RLS		
Adapter, C	losed		
Cap Size	Hose Barb Size	No. of Quick Connects	Qty/pk
53B	n/a	n/a	1

205-533	1-RLS		
Adapter, <sup>1</sup> /	8" HB x 2 and	vent	
Cap Size	Hose Barb Size	No. of Quick Connects	Qty/pk
53B	1/8"	n/a	1

M - M	205-533	2-RLS		
	Adapter, <sup>1</sup> /	4" HB x 2 and	vent	
	Cap Size	Hose Barb Size	No. of Quick Connects	Qty/pk
	53B	1/4"	n/a	1

	205-533	3-RLS		
11 n 11/	Adapter, 5	/ <sub>16</sub> " HB x 2 and	l vent	
	Cap Size	Hose Barb Size	No. of Quick Connects	Qty/pk
CION.	53B	5/16"	n/a	1





# VersaCap® Technology foxxlifesciences.com



#### VersaCap® 83B



205-1003-RLS				
Cap, Clos	sed			
Cap Size	Hose Barb Size	No. of Quick Connects	Qty/pk	
83B	n/a	n/a	2	



205-100	4-RLS			
Cap, Open	with closed ac	dapter		
Cap Size	Hose Barb Size	No. of Quick Connects	Qty/pk	
83B	n/a	n/a	2	



205-100	/-RL5				
Cap, Open top with adapter, 1/4" HB x 2 and Vent					
Cap Size	Hose Barb Size	No. of Quick Connects	Qty/pk		
83B	1/4"	n/a	1		



205-1008-RLS					
Cap, Open top with adapter, ½" HB x 2 and Vent					
Cap Size	Hose Barb Size	No. of Quick Connects	Qty/pk		
83B	1/2"	n/a	1		



205-300	1-RLS		
Adapter, C	losed		
Cap Size	Hose Barb Size	No. of Quick Connects	Qty/pk
83B	n/a	n/a	1



205-3002-RLS				
Adapter, 1/2	" HB x 2			
Cap Size	Hose Barb Size	No. of Quick Connects	Qty/pk	
83B	1/4"	2	1	



205-3003-RLS				
Adapter, 3/4	"HB x 2			
Cap Size	Hose Barb Size	No. of Quick Connects	Qty/pk	
83B	3/8"	2	1	



	205-300	4-RLS			
	Adapter, 1/4	" HB x 3			
1	Cap Size	Hose Barb Size	No. of Quick Connects	Qty/pk	
)	83B	1/4"	3	1	

j		3
	۲	3
1		

205-300	5-RLS		
Adapter, <sup>3</sup> /	″a" HB x 3		
Cap Size	Hose Barb Size	No. of Quick Connects	Qty/pk
83B	3/8"	3	1



205-300	6-RLS			
Adapter, 1/4	" HB x 4			
Cap Size	Hose Barb Size	No. of Quick Connects	Qty/pk	
83B	1/4"	4	1	



205-300	7-RLS			
Adapter, 3/4	"HB x 4			
Cap Size	Hose Barb Size	No. of Quick Connects	Qty/pk	
83B	3/8"	4	1	



	205-500	1-RLS			
	Adapter, 1/4	"HB x 2 and V	ent ent		
	Cap Size	Hose Barb Size	No. of Quick Connects	Qty/pk	
7	83B	1/4"	n/a	1	



205-500	Z-KLO			
Adapter, 1/2" HB x 2 and Vent				
Cap Size	Hose Barb Size	No. of Quick Connects	Qty/pk	
83B	1/2"	n/a	1	



208-300	1-RLS		
Adapter, C	losed Clearvie	w®	
Cap Size	Hose Barb Size	No. of Quick Connects	Qty/pk
83B	n/a	n/a	1



## VersaCap® Technology

foxxlifesciences.com



#### VersaCap® 120mm



205-2003-RLS				
Cap, Clos	sed			
Cap Size	Hose Barb Size	No. of Quick Connects	Qty/pk	
120mm	n/a	n/a	2	



205-400	205-4004-RLS				
Adapter, 1/4	" HB x 3				
Cap Size	Hose Barb Size	No. of Quick Connects	Qty/pk		
120mm	1/4"	3	1		



205-2004-RLS					
Cap, Open	with closed ac	dapter			
Cap Size	Hose Barb Size	No. of Quick Connects	Qty/pk		
120mm	n/a	n/a	2		



205-4005-RLS				
Adapter, 3/8	"HB x 3			
Cap Size	Hose Barb Size	No. of Quick Connects	Qty/pk	
120mm	3/8"	3	1	



205-4001-RLS					
Adapter, C	Adapter, Closed				
Cap Size	Hose Barb Size	No. of Quick Connects	Qty/pk		
120mm	n/a	n/a	1		



205-400	205-4006-RLS				
Adapter, 1/4	" HB x 4				
Cap Size	Hose Barb Size	No. of Quick Connects	Qty/pk		
120mm	1/4"	4	1		



205-4002-RLS						
Adapter, 1/2	Adapter, ¼" HB x 2					
Cap Size	Hose Barb Size	No. of Quick Connects	Qty/pk			
120mm	1/4"	2	1			



205-400	7-RLS			
Adapter, 3/4	"HB x 4			
Cap Size	Hose Barb Size	No. of Quick Connects	Qty/pk	
120mm	3/8"	4	1	



205-400	205-4003-RLS				
Adapter, 3/	"HB x 2				
Cap Size	Hose Barb Size	No. of Quick Connects	Qty/pk		
120mm	3/8"	2	1		



205-400	8-RLS			
Adapter, E	Z Top <sup>®</sup> adapter			
Cap Size	Hose Barb Size	No. of Quick Connects	Qty/pk	
120mm	n/a	n/a	1	

HB = Hose Barb



### VersaBarb® and Spigot Technology

foxxlifesciences.com



#### VersaBarb® and Spigot Technology



# 225-4001-RLS VersaBarb®, Spigot fitting, closed cap Hose Barb Thread Qty/pk Size Diameter n/a 1 1/8" 1



225-4005	-RLS		
<b>Versa</b> Barb <sup>®</sup>	, Spigot fitting,	with 3/4" HB	
Hose Barb Size	Thread Diameter	Qty/pk	
3/4"	1 1/8"	1	



225-4002-RLS					
VersaBarb <sup>®</sup> ,	VersaBarb®, Spigot fitting, with 1/4" HB				
Hose Barb Size	Thread Diameter	Qty/pk			
1/4"	1 1/8"	1			



225-4006	-RLS	
VersaBarb <sup>®</sup> ,	Spigot fitting,	with 3/4" Sanitary Connector
Hose Barb Size	Thread Diameter	Qty/pk
n/a	1 1/8"	1



225-4003	3-RLS	
<b>Versa</b> Barb <sup>®</sup>	Spigot fitting,	with 3/8" HB
Hose Barb Size	Thread Diameter	Qty/pk
3/8"	1 1/8"	1



225-1001-RLS				
Spigot, 2 Pac				
Hose Barb Size	Thread Diameter	Qty/pk		
n/a	1 1/8"	2		



# 225-4004-RLS VersaBarb®, Spigot fitting, with ½" HB Hose Barb Thread Qty/pk Size Diameter ½" 1 ½" 1

HB = Hose Barb

We recommend testing for your specific application before use.



Rev. 1/15/2020 PN: 225-001

