



- Manufactured to accept additional fluid transfer components such as aseptic connectors and air vent filters
- Flexible system is designed for closures that fit bottles or carboys with 53B, 83B, 120mm necks
- Eliminate pouring of high valued products through this leak-proof aseptic system
- Tubing manufactured from high purity silicone
- Polypropylene **VersaCap®**
- All components are autoclavable
- **Made in the USA** 



ISO 9001
ISO 13485


Foxx Life Sciences' **VersaCap®** Filling and Venting closures are aseptic handling solutions for pharmaceutical, biotechnology, and laboratory liquid transfer applications. All closures can be used on bottles and carboys with 53B, 83B, or 120mm closures. The pre-mounted lengths of platinum-cured silicone tubing allow for easy customization of carboys and bottles.




1 Multiple hose barb sizes available

2 **VersaCap®** technology

3 Pre-mounted lengths of platinum-cured silicone





Item Number	Description	Closure Size	Tubing Diameter	Tubing Length	Pack Size
205-6001-RLS	VersaCap® with molded 1/8" hosebarb adapters	53B	1/8"	12"	1
205-6003-RLS	VersaCap® with molded 1/4" hosebarb adapters	53B	1/4"	12"	1
205-6103-RLS	VersaCap® with molded 1/4" hosebarb adapters	83B	1/4"	27"	1
205-6107-RLS	VersaCap® with molded 1/2" hosebarb adapters	83B	1/2"	27"	1
205-6203-RLS	VersaCap® with molded 1/4" hosebarb adapters	120mm	1/4"	37"	1
205-6207-RLS	VersaCap® with molded 1/2" hosebarb adapters	120mm	1/2"	37"	1



BPA -Free

