



We put space to work.

Item # _____

Job _____

MetroBasix™

Reliable performance for any application.

MetroBasix is constructed from a strong steel alloy to endure the rigors of a daily routine while providing a safe, ergonomic and flexible work area for the user.

- **Safe** — The polymer base protects the cart and surrounding environment from incidental contact. The base is designed to allow the 5" (127mm) casters to be mounted outboard of the cabinet, which keeps the cart stable when a fully loaded drawer is accessed.
- **Ergonomic** — Polymer top provides a user-friendly work surface. Optional accessory push handles can be mounted in an ergonomic fashion to allow for easy pushing or pulling of the cart.
- **Flexible** — User defined configurations can be accommodated by filling an empty cart shell with desired 3", 6", 9" and 12" (76, 152, 229, 305 mm) drawer heights.
- **Accessible** — Soft close, full-extension drawers and a full array of accessories put equipment and supplies where needed for convenient access. Optional drawer dividers keep supplies neat and organized.
- **Mobile** — This cabinet comes with 5" (127mm) diameter casters (2 with swivel brakes, 2 with swivel-lock) allowing it to be rolled easily from location to location or moved for convenience when cleaning.
- **Accessories** — A full array of task support accessories are available for above the work surface or the side of the cart.



MBX3210TB-R

(Shown with optional accessories: IV pole, cart top rail, O₂ holder, and push handle.)



MBX3210ANES3

All Metro Catalog Sheets are available on our website: www.metro.com



InterMetro Industries Corporation
 North Washington Street, Wilkes-Barre, PA 18705
 Product Information. U.S. and Canada: 1.800.992.1776
 Outside U.S. and Canada: www.metro.com/contactus

L03-244
Printed in U.S.A. Rev. 1/24

Information and specifications are subject to change without notice. Please confirm at time of order.

Copyright © 2024 InterMetro Industries Corp.




MetroBasix™

36-20



Basix Part Numbering

MBX	Drawers	Body	Lock Options	Color
Metro Basix	Position 1: # of 3" Drawers	S = Solid Color	N = Non-Locking	See Below
	Position 2: # of 6" Drawers	T = Two-Tone	L = Key Lock	
	Position 3: # of 9" Drawers		B = Breakaway Bar	
	Position 4: # of 12" Drawers		P = Push Button Lock	
			E = Electronic Lock	
			R = Prox Lock	



BG
Blue Gray



DB
Dark Blue




TP
Taupe




LB
Light Beige



R
Red



SB
Slate Blue



Y
Yellow

Description	Weight Capacity*		Dimensions		Approx. Pkd. Wt.		Cat. No.
	(lbs.)	(kg)	Height/Width/Length (in.)	(mm)	(lbs.)	(kg)	
3" (76 mm) Deep Steel Drawer	50	24	2.88x17x22.25	73x432x565	7	3.18	MBX103
6" (152 mm) Deep Steel Drawer	50	24	6.25x17x22.25	159x432x565	6	2.73	MBX106
9" (229 mm) Deep Steel Drawer	50	24	9.38x17x22.25	238x432x565	6	2.73	MBX109
12" (318 mm) Deep Steel Drawer	50	24	12.5x17x22.25	318x432x565	6	2.73	MBX112

*Drawer capacity is determined by cart stability with 1 drawer fully extended.

Example Part Numbering

MBX	Drawers	Body	Lock Options	Color	
Metro Basix	3110 Position 1: (x3) 3" Drawers Position 2: (x1) 6" Drawers Position 3: (x1) 9" Drawers Position 4: (x0) 12" Drawers	S Solid	L Key Lock	- Y Yellow	= MBX3110SL-Y



*Accessories shown are separate line items not included in the base cart.



We put space to work.

Item # _____

Job _____



MetroBasix™

Specifications:

Cart Body: Double wall construction
Inner wall: 22 gauge (.030"/.8mm) steel
Outer wall: 20 gauge (.036"/.9mm) steel
Finish: Powder coated

24" (610mm) Body
Height: 38.5" (978mm)
Width: 24.13" (613mm)
Length: 30" (762mm)
Weight: 151 lbs. (72kg)
Drawer Space: 24" (610mm)

30" (762mm) Body
Height: 44" (1140mm)
Width: 24.13" (613mm)
Length: 30" (762mm)
Weight: 180 lbs. (82kg)
Drawer Space: 30" (762mm)

Base: Blow-molded, high-density polyethylene
Steel reinforced caster channels
Caster mounts located outside of cart body footprint to provide counterbalance to full extension drawers

Work Surface: .125" (3mm) ABS , raised on four sides
Dimensions: 1¹/₄" H x 18¹/₂" W x 26³/₄" L (32H x 470W x 679mm L)
Usable work space: 17" x 25" (432mm x 635mm)

Casters: 5" (127mm) diameter with polyurethane tread
2 swivel-brake, 2 swivel-lock (coverts from swivel to rigid)
200 lbs. (90kg) per caster, 600lbs. (270kg) per cart

Drawers: 50 lbs. (23kg) load rating per drawer

Body: 22 Gauge (.030"/.8mm) steel w/16 gauge (.060"/2mm) lock clips
Finish: Powder coated

Fronts: .032" (.8mm) coated aluminum
Finish: Powder coated

Pull: Polymer

Slides: Soft Close
Zinc-Plated Steel
Load Rating: 100 Lbs.(45kg)/set

Lock: Bright Nickel Type 3
G Cam, 7/8" (22mm) Cylinder
2 Keys,
Coded on lock and keys

an Ali Group Company



The Spirit of Excellence

36-20