

## Adapters for Exhaust Filter Connections

### Space problems? Offset adapter!

Space problem in the lab or on small Safety Waste Caps? Not an issue – attach the exhaust air filter to the waste containers in any position. Practical, space-saving and flexible. With the extension **A**, you can create more freedom of movement when attaching tubes. The adapters **A B C D** can be combined with each other.



### Usage with other systems!

You are already using a system e.g. from VICI or Phenomenex? This adapter allows you to use SCAT exhaust filters.



**A** 107 621



**B** 107 624



**C** 107 627



**D** 107 622



Adapter UNF 1/4"-28G to GL 14

Blind plug

**E** 107 620

- Simply closes the exhaust air connection when not in use



**F** 107 632

- Connectivity for HPLC capillaries to the exhaust air connection



Fig.	Part No.	Description	Thread 1	Thread 2	Material
<b>A</b>	<b>107 621</b>	Offset adapter, extension for exhaust filter, GL 14 (f) to GL 14 (m)	GL 14 (f)	GL 14 (m)	PE-HD-EL
<b>B</b>	<b>107 624</b>	Offset adapter 90°, for exhaust filter, GL 14 (f) to GL 14 (m)	GL 14 (f)	GL 14 (m)	PE-HD-EL
<b>C</b>	<b>107 627</b>	Offset adapter 45°, for exhaust filter, GL 14 (f) to GL 14 (m)	GL 14 (f)	GL 14 (m)	PE-HD-EL
<b>D</b>	<b>107 622</b>	Offset adapter 90°, long, for exhaust filter, GL 14 (f) to GL 14 (m)	GL 14 (f)	GL 14 (m)	PE-HD-EL
<b>E</b>	<b>107 620</b>	Blind plug for the exhaust filter connection	GL 14 (m)	-	PTFE
<b>F</b>	<b>107 632</b>	Adapter capillary connection to exhaust filter connection	GL 14 (m)	UNF 1/4" 28G (f)	PTFE
<b>G</b>	<b>108 181</b>	Exhaust filter adapter, suitable for connecting a SCAT exhaust filter on a Waste Cap of third party manufacturers (e.g. Vici)	GL 14 (f)	UNF 1/4" 28G (m)	PE-HD