

## Adapter for JUSTRITE® Containers

### SCAT Europe Waste Systems fit on Containers of JUSTRITE®

Proven SCAT Europe safety at JUSTRITE® containers.  
Suitable adapters and CPC®- couplings.

#### For quick couplings made of plastic



#### For quick couplings made of stainless steel



Fig.	Part No.	Description	Capillary connections	Tube connections	Connection for Exhaust Filter	Material	Unit
	107 628	4-way-collector for plastic coupling	3x	1x	-	PTFE / PFA / PP	1
<b>A</b>	160 527	12-way-collector for plastic coupling	6x	6x	-	PP	1
<b>B</b>	350 099	LISA for plastic coupling	4x	3x	1x	PTFE	1
<b>C</b>	318 957	MARCO for plastic coupling	-	-	-	PE-HD	1
<b>D</b>	107 617	Adapter SCAT Exhaust Filter to plastic coupling	-	-	1x	PE-HD-EL	1
<b>E</b>	160 524	12-way-collector for steel coupling	6x	6x	-	PTFE-EL / stainless steel	1
<b>F</b>	107 631	4-way-collector for steel coupling	3x	1x		PTFE	1
<b>G</b>	107 610	Adapter SCAT Exhaust Filter to steel coupling	-	-	1x	PE-HD-EL	1

HPLC-Supply

HPLC-Disposal

Safety Funnels

Level Control

Containers

Accessory

Addendum