


TUBING AND FITTINGS

Sani-Tech®'s complete line of sanitary tubing and fittings is ideal for high-purity systems. Sanitary flange connections offer the purity required for the most demanding applications.

These products are manufactured with high-purity resins that meet FDA and USP standards and regulations. *Select from the following parts:*



TUBING

Sani-Tech® polypropylene and PVDF tubing is extruded to stainless steel tubing dimensions to allow smooth connection to stainless steel components with no pressure drop.

TEES



Available in standard, short bull, cross tee, reducing, and threaded styles.

ELBOWS



Choose from 45° sweep elbows and 90° sweep and sharp elbows.

SANITARY REDUCERS



Choose from concentric and eccentric styles that connect any size tubing to any smaller size.

END CAPS



Solid and tapped styles available. Your choice of tap size.

SANITARY TRUE UNION CLAMPS



Standard, heavy-duty and ultra heavy-duty styles available.

GASKETS



Choose from silicone, Buna-N, EPDM and Viton®.



All Sani-Tech® fittings are also offered with S.I.B.® ends for use with our S.I.B.® welding equipment.

