

## SUITABLE FOR USE IN CONTROLLED ENVIRONMENTS



### FEATURES & BENEFITS

ComforTech frocks provide an excellent filtration efficiency rate and protect wearers from light chemical splash, particles and biologicals while keeping them cool and comfortable. These frocks are constructed using a high quality microporous composite fabric that is breathable, lightweight and strong. The ComforTech frock is also durable and will effectively meet your controlled environment needs.

A mandarin style snap collar adds style and comfort while providing additional coverage for the chest and neck area. ComforTech frocks are cut and sewn to an appropriate length to ensure full coverage and with the nylon zipper option, which adds another level of particle holdout and protection.

#### Construction Details:

Two styles to choose from: snap or zipper close. Both have elastic wrist with a snap mandarin collar, no pockets and strong, stress-resistant serged seams.

#### Styles/Sizes

Part No.	Size
FK-J2121-1	S
FK-J2121-2	M
FK-J2121-3	L
FK-J2121-4	XL
FK-J2121-5	2X
FK-J2121-6	3X
FK-J2121-7	4X
FK-J2121-8	5X

Raglan sleeve, elastic wrist, snap collar, snap close, serged seams, white, 25/case

#### Styles/Sizes

Part No.	Size
FK-J2141-1	S
FK-J2141-2	M
FK-J2141-3	L
FK-J2141-4	XL
FK-J2141-5	2X
FK-J2141-6	3X
FK-J2141-7	4X

Raglan sleeve, elastic wrist, snap collar, zip close, serged seams, white, 25/case

ComforTech® Frocks	Size	Size Code	Body (A)	Chest (B)	Sleeve (C)
	S	-1	40 ½"	21 ¾"	33 ½"
	M	-2	42 ½"	22 ¾"	34 ½"
	L	-3	43 ½"	23 ¾"	35 ½"
	XL	-4	46"	24 ¾"	36 ½"
	2X	-5	47 ½"	26 ½"	37 ½"
	3X	-6	51 ½"	29"	40"
	4X	-7	52 ½"	29 ½"	41"
	5X	-8	53 ½"	30"	42"