

## SUPERIOR ANTI-SKID PROPERTIES

### FEATURES & BENEFITS

AquaTrak® MaxGrip® Boot Covers provide slip and fall protection on damp and dry floors. With the MaxGrip® anti-skid soles in these top quality boot covers, safety is paramount. The sole material is also extremely clean, durable and strong enough to stand up to the most abrasive floor surfaces. The AquaTrak® upper material is soft, strong, flexible and a fluid resistant material that also provides appropriate particle holdout. These high top boot covers fit over the calf to provide a comfortable fit and prevent the boot cover from sliding down.

#### Construction Details:

Two designs to choose from: elastic top and ankle, and attached ties at ankle for secure fit. Different thread colors around the opening provide convenient size identification.



#### Elastic Ankles

##### Styles/Sizes

Part No.	Size	Boot Color	Thread Color
BT-T5W15-B	M	White	Yellow
BT-T5W12-B	Univ.	White	White
BT-T5W13-B	XL	White	Black

Elastic top and ankle, anti-skid sole, serged seams, 100 pair/case

#### Ankle Ties

##### Styles/Sizes

Part No.	Size	Boot Color	Thread Color
BT-T4W15-B	M	White	Yellow
BT-T4W12-B	Univ.	White	White
BT-T4W13-B	XL	White	Black

Elastic top with ankle ties, anti-skid sole, serged seams, 100 pair/case

MaxGrip® Boot Covers w/Elastic Ankles	Size	Length (A)	Height (B)	Top Relaxed (C)	Top Stretched (C)	Ankle Relaxed (D)	Ankle Stretched (D)	Sole Width (E)
	M	10 ¾"	17"	5 ½"	9"	6"	9 ¾"	4 ¼"
	Univ.	13 ½"	17 ¼"	5 ¾"	9 ½"	6 ½"	11"	4 ½"
	XL	15"	18"	6"	11 ½"	6 ¾"	12 ½"	4 ¾"

MaxGrip® Boot Covers w/Ankle Ties	Size	Length (A)	Height (B)	Top Relaxed (C)	Top Stretched (C)	Tie Length (D)	Sole Width (E)
	M	10 ¾"	18"	4 ½"	9"	17"	4 ¼"
	Univ.	13 ½"	18 ½"	5 ¼"	9 ½"	17"	4 ½"
	XL	15"	19"	5"	11 ½"	17"	4 ¾"