

MULTIPLE STYLES & PACKAGING OPTIONS FOR MAXIMUM VERSATILITY



FEATURES & BENEFITS

Proven in the marketplace, SureGrip shoe covers have established themselves as the shoe cover of choice for a wide range of floor surfaces and in many different applications. Combining durability and cleanliness with outstanding anti-skid properties, SureGrip provides you with a blend of unmatched performance and overall value.

Multiple style and packaging options add to the versatility of SureGrip: serged seams for added strength, traction and a snug fit for safety; unique conductive shoe cover for those environments where electro-static dissipation (ESD) protection is required. The Gaylord packaging option (30 individual bags of 300 each per master carton) reduces waste, labor, and purchase orders while increasing your shelf-space and providing a reusable, pallet sized container for your recyclables or other materials.



Construction Details:

Two styles to choose from: serged seam and conductive to meet any/all market applications. Different thread colors around the opening provide convenient size identification.

Serged Seams Styles/Sizes

Part No.	Size	Shoe Color	Thread Color
SH-X1222-B	Univ.	Blue	White
SH-X1223-B	XL	Blue	Black
SH-X1224-B	2X	Blue	Blue

Non-skid, serged seams, 150 pair/case

Conductive Styles/Sizes

Part No.	Size	Shoe Color	Thread Color
SH-12522-B	Univ.	Blue	White
SH-12523-B	XL	Blue	Black

Conductive for ESD concerns, anti-skid, serged seams, 150 pair/case

Gaylord Packaging Styles/Sizes

Part No.	Size	Shoe Color	Thread Color
SH-X1272-B	Univ.	Blue	White
SH-X1273-B	XL	Blue	Black
SH-X1274-B	2X	Blue	Blue

Non-skid, serged seams, 4,500 pair/case

SureGrip® Serged Seam Shoe Covers	Size	Length (A)	Height (B)	Top Relaxed (C)	Top Stretched (C)
	Univ.	15 ½"	6 ¼"	4 ½"	9 ½"
	XL	17"	6 ½"	5"	11 ½"
	2X	18 ½"	6 ¾"	5 ½"	12 ½"

SureGrip® Conductive Shoe Covers	Size	Length (A)	Height (B)	Top Relaxed (C)	Top Stretched (C)
	Univ.	15 ½"	6 ¼"	4 ½"	9 ½"
	XL	17"	6 ¾"	5"	11 ½"