

Critical Cover® AlphaAir® & AlphaGuard® Face Masks



CRITICAL COVER® ALPHAIR® FACE MASKS *PARTICLE FILTRATION TO 0.1 MICRON*

FEATURES & BENEFITS

The high value AlphaAir® 7" fluid resistant facemasks are clean, comfortable and provide outstanding particle filtration down to 0.1 micron. The full-length aluminum nosepiece helps provide a perfect fit every time. The sonic welded design ensures maximum protection for the user. Multiple colors available and two styles to choose from; soft knit ear loops and ties ultrasonically attached to the mask. ASTM F2100 Compliant Level 3 mask.

AlphaAir Tie Masks Styles/Sizes		AlphaAir Ear Loop Masks Styles/Sizes	
Part No.	Color	Part No.	Color
BL 5100	Blue	BL 5005	Blue
300/case (6 boxes of 50)		GR 5005	Green
		PK 5005	Pink
		YL 5005	Yellow
		500/case (10 boxes of 50)	
AlphaAir Tie Masks Styles/Sizes			
Part No.	Color		
BL 5500	Blue		
500/case (10 boxes of 50)			



CRITICAL COVER® ALPHAGUARD® FACE MASKS *SUITABLE FOR USE IN CONTROLLED ENVIRONMENTS*

FEATURES & BENEFITS

For superior comfort and protection, the AlphaGuard® facemask exhibits an excellent water vapor transmission rate and particle filtration down to 0.1 micron. This 7" fluid resistant facemask has soft, non-irritating knit ear loops that are ultrasonically welded to the mask. ASTM F2100 Compliant Level 3 mask.

AlphaGuard Ear Loop Masks Styles/Sizes	
Part No.	Color
9050	Blue
Cleanroom Packaging 500/case (10 bags of 50)	