

# Critical Cover® CoolOne™ Face Masks



## CRITICAL COVER® COOLONE™ FACE MASKS COOL DESIGN, COOL COMFORT

### FEATURES & BENEFITS

Thanks to the patented Magic Arch® support, the 7" fluid resistant CoolOne face masks give users cooler, more comfortable performance than other designs currently on the market. The integrated Magic Arch® support strip holds the mask away from the nose and mouth to create a larger breathing chamber to enhance user comfort. Multiple colors and three styles to choose from: soft ties which are ultrasonically welded to the mask, comfortable soft knit ear loops or soft knit ear loops ear loops with an anti-fog foam strip. The anti-fog version is for use in high humidity environments or with protective eyewear. Made in the USA. ASTM F2100 Compliant Level 1 mask.

#### CoolOne Tie Masks Styles/Sizes

Part No.	Color
<b>BL 6001</b>	Blue
300/case (6 boxes of 50)	

#### CoolOne Tie Masks Styles/Sizes

Part No.	Color
<b>BL 6001 CC</b>	Blue
Cleanroom Packaging 300/case (6 bags of 50)	

#### CoolOne Ear Loop Masks Styles/Sizes

Part No.	Color
<b>BL 6055</b>	Blue
<b>GR 6055</b>	Green
<b>PK 6055</b>	Pink
<b>YL 6055</b>	Yellow
500/case (10 boxes of 50)	

#### CoolOne Anti-Fog Ear Loop Masks Styles/Sizes

Part No.	Color
<b>BL 6155</b>	Blue
500/case (10 boxes of 50)	



## CRITICAL COVER® COOLONE™ SENSITIVE SKIN FACE MASKS FOR SENSITIVE SKIN

### FEATURES & BENEFITS

If you have skin sensitivity, choose the CoolOne Sensitive Skin face mask. These 7" fluid resistant face masks are free of pigments and dyes. The patented Magic Arch® support strip holds the mask away from the nose and mouth to create a larger breathing chamber to enhance user comfort. Two styles to choose from: soft knit ear loops for a secure and comfortable fit or the anti-fog version for use in high humidity environments or with protective eyewear. Made in the USA. ASTM F2100 Compliant Level 1 mask.

#### CoolOne Sensitive Skin Ear Loop Masks Styles/Sizes

Part No.	Color
<b>WH 1155</b>	White
500/case (10 boxes of 50)	

#### CoolOne Sensitive Skin Ear Anti-Fog Ear Loop Masks Styles/Sizes

Part No.	Color
<b>WH 6155</b>	White
500/case (10 boxes of 50)	