

Critical Cover® ShieldMate® Masks & Shields

CRITICAL COVER® SHIELDMATE® MASK WITH EYE SHIELD



FEATURES & BENEFITS

Get eye and face splash protection and inhalation security in this top value design. The ShieldMate combination facemask and eye shield features a low glare, anti-fog treated, optically clear shield ultrasonically welded to the mask. An anti-fog foam strip on the mask eliminates eyewear fogging. The patented Magic Arch™ support strip holds the mask away from the nose and mouth to create a larger breathing chamber without limiting downward vision. This results in cooler, more comfortable performance for the wearer. ASTM F2100 Compliant Level 1 masks.

ShieldMate Mask with Eye Shield, Ties

Styles/Sizes

Part No.	Color
----------	-------

BL 8101 S	Blue
------------------	-------------

100/case (4 boxes of 25)

ShieldMate Mask with Eye Shield, Ear Loop

Styles/Sizes

Part No.	Color
----------	-------

BL 6155 S	Blue
------------------	-------------

100/case (4 boxes of 25)



CRITICAL COVER® SHIELDMATE® MASK WITH EXTENDED SHIELD

FEATURES & BENEFITS

Get both eye and face splash protection in this truly unique design. This ShieldMate® is a 100% polypropylene mask, offering .1 micron filtration. The mask design features comfortable ear loops and a full-length aluminum nosepiece. The 5-mil shield is low glare, optically clear and distortion-free. The wraparound shield has permanent anti-fog treatment on both sides. ASTM F2100 Compliant Level 3 mask.

ShieldMate Mask with Extended Shield, Ear Loop

Styles/Sizes

Part No.	Color
----------	-------

SM 2005X	Blue
-----------------	-------------

100/case (4 boxes of 25)