

## Lab-Tips® - LTN1465 Long Handle Nonwoven Polyester Swab

100% hydroentangled polyester long handle swab

The LTN1465 long handled nonwoven polyester swab is a general purpose swab specially suited for cleaning hard to reach areas. Its flexible paddle tip and soft non-abrasive surface make it a good choice for sensitive applications.



### Other Nonwoven tip swabs

- LTN70F Small handle

### Key Attributes

- Double layer, hydro-entangled polyester tip over flexible paddle
- Soft and non abrasive
- Thermally bonded head
- Polypropylene handle
- Packaged in an ISO class 4 (Class 10) cleanroom

### Benefits

- Low particle generation
- Low non volatile residue
- Excellent chemical resistance
- No contaminating adhesives
- Good sorbency
- Free from Silicone, amides and DOP

### Environmental

- Complies with RoHS and REACH requirements

### Applications

- Cleaning of sensitive equipment
- Applying lubricants and other liquids
- Cleaning validation
- Surface sampling
- Cleaning with solvents

### Physical Characteristics

<b>Head Material</b>	Nonwoven Polyester	<b>Handle Thickness</b>	3.1mm {0.157"}
<b>Head Width</b>	7mm {0.275"}	<b>Handle Length</b>	145mm {5.71"}
<b>Head Thickness</b>	4mm {0.118"}	<b>Total Swab Length</b>	162mm {6.38"}
<b>Head Length</b>	17mm {0.670"}	<b>Head Bond</b>	Thermal
<b>Handle Material</b>	Polypropylene	<b>Handle Text</b>	Berkshire
<b>Handle Width</b>	3.1mm {0.157"}		

[www.berkshire.com](http://www.berkshire.com)

Contact: Tel 1 800 242 7000

[info@berkshire.com](mailto:info@berkshire.com)

Americas	Tel 1 413 528 2602	Fax 1 413 528 2614	<a href="mailto:info@berkshire.com">info@berkshire.com</a>
Europe	Tel 44 (0) 870 757 2877	Fax 44 (0) 870 757 2878	<a href="mailto:enquiries@berkshire.uk.com">enquiries@berkshire.uk.com</a>
SE Asia	Tel 65 6252 4313	Fax 65 6252 4312	<a href="mailto:enquiries@berkshire.com.sg">enquiries@berkshire.com.sg</a>
Japan	Tel 81 3 5827 2380	Fax 81 3 5827 2382	<a href="mailto:master@berkshire.co.jp">master@berkshire.co.jp</a>

## Technical Data:

ATTRIBUTE		LTN1465	TEST METHOD
Non-Volatile Residue	IPA mg/swab	0.060	IEST-RP.CC004.3, Sec 7.1.2
Ions (µg/swab)	Chloride	0.10	IEST-RP.CC004.3, Sec 7.2.2
	Sulfate	0.077	
	Nitrate	ND*	
	Phosphate	ND*	
	Fluoride	ND*	
	Potassium	0.0088	
	Calcium	0.028	
	Sodium	0.027	
	Magnesium	0.0057	

## Notes:

- Data shown are typical values
- Due to differences in test methods applied and equipment utilized by different swab manufacturers, valid product comparisons may only be obtained through side-by-side testing in the same test facility, under similar conditions.
- \* ND\* = below the detection limit of the equipment.

## Order Information:

Product	Part Number	Size	Swabs/Pk	Pks/Cs	Case size
Long Handle Nonwoven	LTN1465.10	162mm	100	10	180 x 150 x 170mm

## Other Berkshire Products



Wipers




Mop Systems



Face Masks



Glove Liners



Documentation Systems



Swabs