



# Transport materials safely & efficiently.

## Introducing Metro's NEW Carboy Cart

Elevate your laboratory operations with Metro's Carboy Cart. Engineered to meet the high standards of modern lab environments by leveraging the advanced MetroMax i design. This cart combines durability, ease of maintenance, and superior mobility to streamline your workflow. Ideal for handling and transporting carboys, it offers a single-shelf design that maximizes space efficiency and facilitates effortless organization.

- MetroMax i shelf and polymer post have a lifetime warranty against corrosion.
- Ergonomic handle
- Weight capacity 750lbs.
- Total Guard all-polymer and stainless casters
- Size: 24"x 30"



\*Shown with 2 50L containers



Holds up to **750** lbs.

# Lab Utility Carts

## High Performance, Corrosion Protection

Durable materials and total-guard polymer casters for long life in laboratory environments.



Solid Shelf

### Metro Carboy Cart

Model No.	Nominal Size		Actual Outside Dimensions (incl. bumpers)		# of Tiers	Approx. Pkd Weight	
	(in.)	(mm)	(in.)	(mm)		(lbs.)	(kg)
<b>CRLS421CB</b>	24x30	457x760	26.3125 x 34.5	668 x 876	1	46.6	21.1

Weight Capacity : 750lbs.

### MetroMax i all polymer carts

- MetroMax i shelves and polymer posts have a lifetime warranty against corrosion.
- Ergonomic Extended Handle
- 2-shelf and 3-shelf models
- Total-guard all-polymer and stainless casters (2 swivel and 2 total lock)
- Choice of shelf styles –
  - Models with all solid shelves
  - Models with a top solid shelf and lower open grid shelves.



Solid Shelves

### MetroMax i Models with All Solid Shelves

Model No.	Nominal Size		Actual Outside Dimensions (incl. bumpers)		# of Tiers	Approx. Pkd Weight	
	(in.)	(mm)	(in.)	(mm)		(lbs.)	(kg)
<b>MXUC1830F-25L</b>	18x30	457x760	20.3125 x 34.5	516 x 876	2	45.7	20.7
<b>MXUC2436F-25L</b>	24x36	610x914	26.3125 x 40.5	668 x 876	2	54.5	24.7
<b>MXUC1830F-35L</b>	18x30	457x760	20.3125 x 34.5	516 x 876	3	58.7	26.6
<b>MXUC2436F-35L</b>	24x36	610x914	26.3125 x 40.5	668 x 876	3	72.5	32.9

Overall Height: 39.25" H (997mm)

### MetroMax i Models with Top Solid Shelf. (Lower tiers are open grid)

Model No.	Nominal Size		Actual Outside Dimensions (incl. bumpers)		# of Tiers	Approx. Pkd Weight	
	(in.)	(mm)	(in.)	(mm)		(lbs.)	(kg)
<b>MXUC1830FT-25L</b>	18x30	457x760	20.3125 x 34.5	516 x 876	2	43.7	19.8
<b>MXUC2436FT-25L</b>	24x36	610x914	26.3125 x 40.5	668 x 876	2	51.5	23.4
<b>MXUC1830FT-35L</b>	18x30	457x760	20.3125 x 34.5	516 x 876	3	54.7	24.8
<b>MXUC2436FT-35L</b>	24x36	610x914	26.3125 x 40.5	668 x 876	3	66.5	30.2

Overall Height: 39.25" H (997mm)



Solid & Grid Shelves

## Polymer/Stainless Steel Total-Guard Stem Casters

Designed to provide excellent corrosion protection and moisture resistance for laboratory and clinical environments.

Non-marking thermoplastic elastomer (TPE) tread (70 shore) is quieter and more shock absorbent than a traditional polyurethane (TPU) tread.

Swivel and total lock brake versions are available.



We put space to work.