

LABORATORY SPILL TRAYS



Compatible with Metro Shelving Products



LABORATORY SPILL TRAYS

- Protects cabinets and shelving from acid spills
- Over 250 sizes of trays available
- Available in 1", 1-1/2" & 2" depths
- Standard color is natural, special colors available
- 3/16" tray walls for rigidity
- One piece rotationally molded polyethylene
- Rounded corners for easy cleaning
- Dishwasher safe
- Compatible with Metro products

TECHNICAL INFORMATION

Laboratory spill trays are a minimum 3/16" thick, one piece stress free, virgin polyethylene with smooth interior walls and covered corners for easy cleaning. Tray conforms to ASTM F 1412 for corrosion resistance, and will provide secondary containment to comply with OSHA 1910.1450 and EPA 264.175 for chemical stockrooms/storerooms.

