


Datasheet

Part No.:	106484
Description	SafetyWasteCap S90, electronic level control, thread GL25 (m), PE-HD electrostatic conductive, ATEX-compliant + connector for exhaust filter
Measurements (WxDxH)	104 x 104 x 180 mm
Weight:	0,657 kg
Colour	black
Material:	PEHD-EL
Sealing material:	PEHD-EL
Chemical resistance:	good
Surface resistance:	$\leq 10^6$ Ohm
Operating temperature:	$\leq 60^\circ$ C
Closure type:	screw cap
Thread size:	S 90

Compliance / approvals		
Document No.:	Date	Inspection authority:
10ATEX 555832 X	11.01.2012	TÜV Nord

UN approval:

Other approvals / marks:
 II 2G Ex ia IIB T6 Gb <ul style="list-style-type: none"> waste: S90 (f) to GL25 (m) vapor: GL14 (f) electronic: M8x1 3-pin (m)

