

A
306 509
4-in-1 Collector
NICOLE
• 4x GL45
• L1 = GL45




HPLC-Multiple collector

The 4-in-1 collector named NICOLE allows you to place up to 4 Safety Waste Caps on a collected drain.

- PE-HD electrostatic conductive
- Can be combined with Safety Waste Caps, table and wall feedthroughs

B
106 584
Thread adapter
• L1 = GL25 (m)
• L2 = G 1 1/2" (w)



Thread adapter

Use the thread adapter to connect the lab sink to the SymLine® disposal system.

- PTFE electrostatic conductive
- Can be combined with pipe system, tube system with SymLine® FlexTube and FlexTube Pro

Fig.	Part No.	Description	Material
A	306 509	4-in-1 collector, NICOLE, 4x GL45 (m) to 1x GL45 (f), incl. blind cap for HPLC connection	PE-HD-EL
B	106 584	Thread adapter for lab sink, PTFE electrostatic conductive, GL25 (m) to G 1 1/2" (f)	PTFE-EL
	106 434	Drain adapter for lab sink, G 1 1/2" (f) / 32 mm OD	PTFE-EL