


## Capillary Connection - Safety Caps - Preparative HPLC Fittings, Blind Plugs

<p><b>A</b> 107 045</p> <p>For Safety Cap:</p> <ul style="list-style-type: none"> <li>■ 308 032</li> <li>■ 309 032</li> <li>■ 310 032</li> </ul>  <p>■ Ø 4.76 mm (UNF 5/16")</p>	<p><b>B</b> 160 503</p> <p>For Safety Cap:</p> <ul style="list-style-type: none"> <li>■ 308 032</li> <li>■ 309 032</li> <li>■ 310 032</li> </ul>  <p>■ Blind plug (UNF 5/16")</p>	<p><b>C</b> 107 044</p> <p>For Safety Cap:</p> <ul style="list-style-type: none"> <li>■ 307 007</li> <li>■ 307 008</li> <li>■ 307 009</li> </ul>  <p>■ Ø 6.35 mm (NPT 1/8")</p>	<p><b>D</b> 160 506</p> <p>For Safety Cap:</p> <ul style="list-style-type: none"> <li>■ 307 007</li> <li>■ 307 008</li> <li>■ 307 009</li> </ul>  <p>■ Blind plug (NPT 1/8")</p>
---	--	---	---



**Good to know!**

With the **SCAT measuring template**, you can easily determine the outside diameter of capillaries and the inside diameter of tubes.

This tool helps you to find the right fittings and tube connectors for your system.

See **page 169** for ordering information.

Fig.	Part No.	Description	Capillary Size Ø OD	Thread	Material	Colour	Unit
	107 047	Fitting for preparative HPLC	4.0 mm	UNF 5/16"	PTFE	White	1
<b>A</b>	107 045	Fitting for preparative HPLC	4.76 mm	UNF 5/16"	PTFE	White	1
<b>B</b>	160 503	Blind Plug for preparative HPLC	-	UNF 5/16"	PTFE	White	10
	160 515	Blind Plug for preparative HPLC	-	UNF 5/16"	PTFE	White	5
	107 046	Fitting for preparative HPLC	6.0 mm	NPT 1/8"	PTFE	White	1
<b>C</b>	107 044	Fitting for preparative HPLC	6.35 mm	NPT 1/8"	PTFE	White	1
<b>D</b>	160 506	Blind Plug for tube connection	-	NPT 1/8"	PTFE	White	1

- HPLC-Supply
- HPLC-Disposal
- Safety Funnels
- Level Control
- Containers
- Accessory
- Addendum